



Instruction for Use

【Product Name】 **UNDERPAD**

【Model】 **Customized**

【 Device description 】 A flat piece of absorbent fabric material [e.g., non-woven polyethylene (PE), cellulose] intended to be placed on a surface (e.g., bed mattress, table, floor) and underneath a patient to absorb bodily fluids resulting from surgical procedures, childbirth and/or incontinence; it does not include an antimicrobial agent(s). It is intended to be used in the healthcare facility and in the home. This is a single-use device.

【 Intended Use 】 The underpads are intended to be placed on a bed mattress surface and underneath a patient to absorb bodily fluids. The product is for single use only, and provided non-sterile.

【 Appearance and Components 】 Made from nonwovens and plastic film composite or sewn. Provide non-sterile, disposable use.

【Usage】

1. Read the label and instruction for use carefully before used;
2. Take off the package, do not use the product once the package was damaged or the product has expired;
3. Put the under pads on the sickbed.

【Precaution and Warning】

- 1.This product is disposable, please timely change the clean product.
- 2.This product is not sterilized before delivery and is not provided aseptically.
3. Do not use it with electric blankets and other heating appliances.

4. This product is not washable. If you wash this product with other clothes, please dehydrate the clothes. After this product residual material off, please carefully wipe the inside of the washing machine after rinsing with water.

5. Please keep the product out of reach of children and old people to avoid accidental ingestion.

【Contraindications】

Do not use if you are allergic to non-woven fabrics



【Storage】

These products should be stored in a cool dry area, away from heat and direct sunlight.

【Shelf Life】

Three years

【Labels, Packing Logo Design】

Symbol	Introductions	Symbol	Introductions
	CE Symbol		medical device
	Manufacturer Name Address		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	Manufacture Date		Use until year & month (Expiration date)
	Batch Code		Do not reuse" are "single use, "Use only once
	Warnings and Precautions		non-sterile

Tianjin Great Hygiene Products Co., Ltd.

Add:Room 1602, Build NO. 1-B, No. 15,Rongyuan Road, Huayuan Industrial Zone,
Binhai High-tech Zone,Tianjin,China

Website: <https://www.greathygiene.com> Whatsapp: +86 15620554861

	Consult instructions for use		Do not use if package is damaged
	Country of Manufacture		Model Number
	Importer		Distributor
	Unique Device Identification		

【Manufacturer】



Company: Tianjin Great Hygiene Products Co. Ltd.
Address:Room 1602, Build NO. 1-B, No. 15,Rongyuan Road,
Huayuan Industrial Zone, Binhai High-tech Zone,Tianjin,China
Tel: 86-13752029469
Email:aaron.zhang@kowford.com

【European Authorized Representative】



HUMISS TRADING S. L.

Calle Luis Buñuel 12, Madrid, 28018, Spain

